

RAID-CRC Screen qPCR Kit

Colorectal Cancer Screening qPCR Kit

INSTRUCTIONS FOR USE

Reference number: **REF**

RAID-CRC Screen qPCR Kit: CRC-01-1230-01

The **RAID-CRC Screen qPCR kit** is an *in vitro* diagnostic device intended exclusively for use by trained laboratory professionals (professional user).

Intended purpose

The **RAID-CRC Screen qPCR kit is intended for the screening of advanced colorectal neoplasia**, including precancerous lesions and colorectal cancer, in asymptomatic individuals aged 50 to 69 years (inclusive) who have obtained a positive result from the faecal immunochemical test (FIT).

The clinical performance of RAID-CRC Screen was assessed in three independent studies conducted in a FIT-positive screening population (Malagón *et al.*, 2019; Malagón *et al.*, 2020). Weighted clinical sensitivity and specificity were 97.2% and 22.1% for colorectal cancer, and 86.8% and 22.1% for advanced colorectal neoplasia. The predictive values were 13.5% (PPV) and 98.3% (NPV) for colorectal cancer, and 48.6% (PPV) and 65.0% (NPV) for advanced neoplasia.

The RAID-CRC Screen assay detects specific bacterial markers in DNA extracted from patient stool samples by means of quantitative polymerase chain reaction (qPCR). The test analyses a panel of **six faecal bacterial markers** indicative of both favourable and unfavourable intestinal health conditions: B46, B48, *Faecalibacterium prausnitzii*, *Gemella morbillorum*, *Bacteroides fragilis*, and Eubacteria.

The kit **enables amplification and quantification of the gene fragments characteristic** of the microorganisms listed above. Results are provided both qualitatively (Positive/Not compatible) and quantitatively (Cq values and genomic copies/ μ L). Result generation is performed through the GoodGut-Test™ platform (<https://goodgut-test.eu>), which has been verified and validated for use with this kit as an *in vitro* diagnostic accessory.

This product is not automated and is intended for use solely by professional laboratory personnel.

Principles of the test

The **RAID-CRC Screen qPCR kit** has been optimised through the analysis of **multiplex qPCR reactions** using specific primers and fluorescent probes. The result is an **easy-to-use** tool with **high sensitivity, specificity, reproducibility**, and a **broad dynamic range**.

The assay is based on the 5'-exonuclease activity of the DNA polymerase enzyme. During DNA amplification, the enzyme cleaves probes hybridised to the complementary DNA sequence, separating the quencher from the reporter. This reaction generates an increase in fluorescence proportional to the amount of hydrolysed DNA sequence. Fluorescence is measured using validated real-time PCR thermocyclers (see Annex 2).

To obtain a complete diagnostic result, each sample requires **three independent qPCR amplifications**: the first, denominated **Screen_A** or **SCR_A**, is multiplex and targets three bacterial markers; the second, **Screen_B** or **SCR_B**, is also multiplex and targets two markers; and the third, **Screen_C** or **SCR_C**, is singleplex and targets a single bacterial marker. **Together, these three amplifications constitute one complete RAID-CRC Screen analytical reaction, corresponding to the full analysis of one sample, whether it is a clinical specimen, a positive control, or a negative control.**



Each kit allows for **230 complete analytical reactions**, with each reaction comprising the Screen_A, Screen_B, and Screen_C amplifications. The **4X Master Mix is provided in 3 vials**, sufficient for 230 complete reactions (approximately 77 reactions per vial). It is supplied ready to use and contains all components required to perform qPCR. **The oligonucleotides (primers and probes) are provided lyophilised in separate vials**, each sufficient for 230 complete reactions. **RNase-free water is supplied in 2 vials**, each allowing for at least 115 complete reactions. Additionally, the kit includes **3 positive controls**, one for each amplification (Screen_A, Screen_B, and Screen_C), in separate tubes, enabling verification of correct qPCR performance.

Note: Throughout this document, the term **reagents** refers collectively to the 4X Master Mix, the oligonucleotides, the RNase-free water, and the positive controls.

Requirements for the RAID-CRC Screen analysis

The **RAID-CRC Screen qPCR kit** has been optimised for the quantification of the bacterial markers specified above from DNA extracted from stool samples. To ensure the validity of the analysis, the following requirements must be met:

Intended test population:

- **Asymptomatic individuals aged 50 to 69 years** (inclusive), considered an **intermediate-risk population** for colorectal cancer.
- The test must be applied exclusively to individuals with a **positive FIT result**, using a **cut-off value of 100 ng of haemoglobin per mL** obtained with the **Eiken Chemical FIT collection tube**, equivalent to 20 µg of haemoglobin per g of faeces.

Note: If a different FIT cut-off is used, please contact the manufacturer (support@goodgut.eu).

Exclusion criteria:

Stool samples must not be collected from individuals who:

- Are pregnant.
- Are undergoing or have received oral or intravenous antibiotic treatment during the month prior to sample collection.
- Have undergone a colonoscopy during the month before sample collection.
- Have previously undergone surgical resection of any part of the digestive tract.
- Have received chemotherapy and/or radiotherapy within the last 6 months.

Pre-analytical requirements:

- FIT determination must be performed using a dedicated analyser from **Eiken Chemical**, and FIT values must be expressed in **ng of haemoglobin per mL**.
- Stool samples **must be processed within 48 hours of collection**.
- During this period, sample collection must be performed using the **Eiken Chemical FIT collection tube**.
- If collection with the FIT tube cannot be performed within the first 48 hours, the sample must be frozen until collection can take place.
- Once collected in the FIT tube, the sample must be stored at 2–8 °C for a minimum of 48 hours prior to FIT determination and subsequent DNA extraction.
- DNA extraction must be performed within 18 days of sample collection using the FIT tube. If this is not possible, the FIT tube must be stored at –20 °C until the day of analysis, after having remained at least 48 hours at 2–8 °C.
- Because the stool sample is diluted in the FIT tube solution, a preliminary preparation step must be performed before DNA extraction:
 1. Homogenise the FIT tube by gentle manual inversion.
 2. Transfer the contents of the FIT tube into a 1.5-mL microtube (typically between 1.0 and 1.5 mL).
 3. Centrifuge the 1.5-mL tubes for 10 minutes at 4,000 × g.

4. Remove the supernatant, leaving 100–200 µL of the initial volume.
5. Homogenise the pellet by pipetting.
6. Transfer the pellet into the bead-containing tubes provided in the DNeasy PowerSoil Pro Kit and proceed according to the manufacturer's instructions.

Note: *Instead of adding 250 mg of soil in step 1 of the PowerSoil protocol, insert the resuspended volume from the FIT tube.*

- After DNA extraction, the DNA may be used immediately for qPCR or stored at –20 °C until analysis.

Compliance with these requirements is essential to ensure that the results obtained fall within the reference ranges established for colorectal cancer screening using the RAID-CRC Screen qPCR kit.

Kit contents

The **RAID-CRC Screen qPCR kit** contains 3 vials of 4X Master Mix, 2 vials of RNase-free water, 18 vials of oligonucleotides (including forward primer, reverse primer, and probe), 3 vials of positive controls, and one Quick Start Protocol leaflet.

Table 1. Components included in the RAID-CRC Screen qPCR kit and reagent information.

RAID-CRC Screen Kit (230 reactions 10 µL/reaction)	Concentration	Vial colour	Quantity
Multiplex Master Mix	4X	Clear vial with red cap	3 vials (625 µL)
Primer SCR_A_f1	2,5 µM	Amber glass	1 vial (80 µL*)
Primer SCR_A_r1	2,5 µM	Amber glass	1 vial (80 µL*)
Probe SCR_A_FAM (containing fluorophore FAM and quencher BHQ1)	2,5 µM	Amber glass	1 vial (120 µL*)
Primer SCR_A_f2	5,0 µM	Amber glass	1 vial (180 µL*)
Primer SCR_A_r2	5,0 µM	Amber glass	1 vial (180 µL*)
Probe SCR_A_HEX (containing fluorophore HEX and quencher BHQ1)	5,0 µM	Amber glass	1 vial (80 µL*)
Primer SCR_A_f3	5,0 µM	Amber glass	1 vial (180 µL*)
Primer SCR_A_r3	5,0 µM	Amber glass	1 vial (180 µL*)
Probe SCR_A_ROX (containing fluorophore ROX and quencher BHQ2)	5,0 µM	Amber glass	1 vial (100 µL*)
Primer SCR_B_f1	5,0 µM	Amber glass	1 vial (180 µL*)
Primer SCR_B_r1	5,0 µM	Amber glass	1 vial (180 µL*)
Probe SCR_B_HEX (containing fluorophore HEX and quencher BHQ1)	5,0 µM	Amber glass	1 vial (80 µL*)
Primer SCR_B_f2	5,0 µM	Amber glass	1 vial (160 µL*)
Primer SCR_B_r2	5,0 µM	Amber glass	1 vial (160 µL*)
Probe SCR_B_CY5 (containing fluorophore CY5 and quencher BHQ2)	5,0 µM	Amber glass	1 vial (120 µL*)
Primer SCR_C_f1	5,0 µM	Amber glass	1 vial (180 µL*)
Primer SCR_C_r1	5,0 µM	Amber glass	1 vial (180 µL*)
Probe SCR_C_FAM (containing fluorophore FAM and quencher BHQ1)	5,0 µM	Amber glass	1 vial (160 µL*)
Positive Control A	10 ⁴ -10 ⁷ copies/µL**	Clear	1 vial (185 µL)
Positive Control B	10 ⁵ -10 ⁶ copies/µL**	Clear	1 vial (185 µL)
Positive Control C	10 ⁸ copies/µL	Clear	1 vial (185 µL)
RNase-free water	NA	Clear	2 vials (1,9 mL)
Quick Start Protocol	NA	NA	1 leaflet

NA: Not applicable. * Volume indicated for resuspending the lyophilised oligonucleotide in Tris-HCl pH 8.0 to obtain a final concentration of 2.5 µM for the multiplex SCR_A marker 1 tubes and 5.0 µM for the remaining markers. ** Depending on the marker.

Reagents, materials and equipment not provided in the kit

To perform the RAID-CRC Screen qPCR analysis correctly, the following reagents, materials, and equipment are required but **not included** in the kit:

- Eiken Chemical FIT Collection Tube
- DNA extraction kit (see Annex 1 for compatibility verification)
- Thermocycler (see Annex 2 for compatibility verification)
- Tris-HCl pH 8.0 buffer (for oligonucleotide resuspension)
- Microcentrifuge tubes
- PCR or qPCR tube strips and optical strip caps (8-tube strips)
- Filtered pipette tips
- Biosafety cabinet
- Mechanical disruptor or vortex with tube adaptor
- Microcentrifuge for 1.5-mL tubes
- Strip-tube centrifuge
- Micropipettes (0.5–10 µL, 10–100 µL, and 100–1000 µL)
- Powder-free disposable gloves

Shipping and storage conditions

The **RAID-CRC Screen qPCR kit is shipped refrigerated** at temperatures between **2 °C and 8 °C**. Upon receipt, the Master Mix and positive controls must be stored between **–30 °C and –15 °C** in a stable-temperature freezer and protected from light. Lyophilised primers and probes may be stored at **room temperature** until they are resuspended in Tris-HCl buffer (pH 8.0). **Once resuspended**, they must be stored between **–30 °C and –15 °C** in a stable-temperature freezer and protected from light.

In-use stability

Shelf-life after opening: once opened, reagents **remain stable until the expiration date indicated on the label**, provided they are stored between **–30 °C and –15 °C**. Outside this temperature range, the product may remain for a **maximum of 24 hours between 2 °C and 8 °C** without affecting its specifications.

Freeze–thaw cycles: reagents may undergo up to **10 freeze–thaw cycles**. If more thawing cycles are expected, it is recommended to prepare additional aliquots of the reagents.

Safety information

- This product is intended **exclusively for professional use**.
- **Do not use** the kit after the expiration date indicated.
- **A unidirectional workflow** must be established, beginning in the extraction area and proceeding to the amplification and detection areas. Samples, reagents, and equipment must **not be returned** to the previous area once each step is completed.
- Prepare the Master Mix in a **DNA-free biosafety cabinet**. **Avoid loading samples and positive controls into the qPCR plate inside the same cabinet**.
- **Good Laboratory Practices (GLP)** must be strictly followed:
 - Wear protective clothing, disposable gloves, safety goggles, and a mask.
 - Do not eat, drink, or smoke in the work area.
 - Wash hands after completing the analysis.
- Consumables and reagents used for qPCR must be disposed of in **biohazard waste containers**.
- **Regular decontamination of equipment**, particularly micropipettes and work surfaces, is recommended.



Caution: The 4X Master Mix contains **1,2,4-triazole**, a hazardous substance. It may impair fertility or cause harm to the unborn child. It may also pose risks to breastfed children. Special instructions must be obtained before use. Do not handle the kit until all safety instructions have been read and understood. Good Laboratory Practices described in this section must be followed, and personnel must be informed of the risks associated with handling the product. **In the event of actual or suspected exposure: consult a physician.**



Caution: Do **NOT** add bleach or acidic solutions directly to sample preparation waste, as toxic irritating vapours may be produced.

Information on interfering substances

Refer to the section **Requirements for the RAID-CRC Screen Analysis** (page 1).

Quality control

In accordance with the **GoodGut Quality Management System** (ISO 13485 certified), each lot of the **RAID-CRC Screen qPCR kit** is tested under predefined specifications to ensure **activity**, **efficiency**, and **sensitivity**.

The Certificate of Analysis can be found in the professional area of the GoodGut website: <https://professionalarea.goodgut.eu/>.

Limitations of use

The reagents in this kit are designed to function exclusively with this qPCR kit. **Use with other assays is not recommended.**

The reagents are only compatible with the instruments specified in **Annex 2**.

To date, no additional components have been identified that could influence the measurements.

Accessories of the RAID-CRC Screen qPCR Kit

To obtain the diagnostic result, use the **GoodGut-Test™** web platform (<https://goodgut-test.eu>). Access to the platform is provided separately upon purchase of the RAID-CRC Screen qPCR kit.

A **user manual** is supplied together with a **DEMO** explaining the operation of the platform, intended for use by professional laboratory personnel. If you have not received it, please contact support@goodgut.eu.

The **recommended computer configuration** for using the GoodGut-Test™ web platform is shown in Table 2.

Table 2. Recommended computer configuration for the GoodGut-Test™ web platform.

	For WINDOWS	For MAC
Scale	125%	125%
Screen resolution	1920 x 1080	1920 x 1080
Screen resolution	Landscape	Landscape

An Internet connection is required to use the GoodGut-Test™ web platform. The platform is compatible with Google Chrome, Microsoft Edge, and Mozilla Firefox.

Reference measurement procedure

To ensure correct performance of the RAID-CRC Screen qPCR kit, **positive controls** of known concentration must be included in every multiplex qPCR run (see *RAID-CRC Screen qPCR Kit Protocol*, page 8). A **No Template Control (NTC)** must also be included to verify the absence of reaction contamination.

Positive control

As the assay is based on a multiplex qPCR that simultaneously analyses 2 or 3 biomarkers, each positive control consists of a mixture of specific target sequences corresponding to each assay. For each lot of positive control, a **tolerance range** has been established based on the analysis of three independent duplicate runs selected at random.

After running and interpreting the assay, **the Ct value of the positive control must fall within the established range for the corresponding lot**. If the Ct value falls outside the accepted range, the results cannot be considered reliable. The **GoodGut-Test™** web platform automatically indicates whether the positive controls are accepted or rejected. In case of rejection, the sample analysis must be repeated.

The **tolerance ranges for positive controls** are available in the technical specifications of the RAID-CRC Screen qPCR kit, supplied with each lot at the time of purchase and accessible in the GoodGut professional area: <https://professionalarea.goodgut.eu/>.

No template control (NTC)

The NTC is used to confirm that the reaction mix is free from contamination. Each marker included in the RAID-CRC Screen qPCR kit has a **predefined minimum Ct value**.

After running the assay, **the Ct value of the NTC must be above the predefined minimum threshold**. If the Ct value of the NTC is below the acceptable minimum, the results are not reliable. The **GoodGut-Test™** web platform indicates whether the NTCs are accepted or rejected. In the event of rejection, the analysis must be repeated.

The minimum accepted Ct values for NTCs are included in the technical specifications of the RAID-CRC Screen qPCR kit, supplied with each lot and accessible in the GoodGut professional area: <https://professionalarea.goodgut.eu/>.

Reagent information

Table 3. Information on the reagents included in the RAID-CRC Screen qPCR Kit.

Component	Description
Multiplex Master Mix 4X	The DNA polymerase is a modified form of a recombinant 94 kDa DNA polymerase isolated from <i>Thermus aquaticus</i> . The enzyme is supplied in an inactive state with no enzymatic activity at room temperature and is activated by incubating for 1 minute at 95 °C. The master mix contains Tris-HCl, KCl, NH ₄ Cl, MgCl ₂ , and additives enabling rapid cycling, as well as ultrapure dATP, dCTP, dGTP, and dTTP.
RNase-free water	Distilled RNase-free water for molecular biology applications.
Primers (forward and reverse)	For each marker, two primers (forward and reverse) are provided, purified by desalting and preloaded in the corresponding vial.
Probes	For each marker, a probe purified by HPLC is provided, preloaded in the corresponding vial.
Positive Control SCR_A Positive Control SCR_B Positive Control SCR_C	Each positive control contains a distinct mixture of qPCR amplification products depending on the specific qPCR assay. All controls undergo rigorous quality verification, including size confirmation by capillary electrophoresis and sequence identification by mass spectrometry.

RAID-CRC Screen qPCR kit protocol

To obtain the RAID-CRC Screen results, the following protocol must be followed.

The concentration of primers/probes, as well as the cycling parameters (temperature for hybridisation, number of cycles, and incubation times), have been fully optimised to achieve optimal analytical performance and specificity.

Before beginning, primers and probes must be resuspended in the volume of Tris-HCl buffer pH 8.0 indicated in the protocol (see Kit Contents section).

Note: For optimal resuspension of primers and probes, after adding the Tris-HCl buffer, incubate the tubes at room temperature for 1 hour or overnight at 4 °C without changing the container. Once resuspended, they must be stored between –30 °C and –15 °C in a stable-temperature freezer and protected from light.

For diagnostic interpretation, two multiplex qPCR analyses and one singleplex qPCR analysis must be performed per sample: **Screen_A (SCR_A)**, **Screen_B (SCR_B)** and **Screen_C (SCR_C)**.

qPCR protocol steps:

1. Preparation of PCR tube strips:

Determine the number of tubes and caps required for the intended reactions, taking into account all samples and controls to be included in each qPCR run (materials not included in the kit). Each qPCR analysis must include one **positive control** and one **No Template Control (NTC)**.

Note: Each qPCR has its own specific positive control.

2. Initial preparation of reagents:

Thaw the 4X Multiplex Master Mix, primers, probes, and positive controls supplied in the qPCR kit.

3. Multiplex qPCR Analysis – Screen A:

Add the components listed in Table 4 into a microcentrifuge tube. It is recommended to prepare a reaction mix volume of $n \times 1.1$, where n is the number of reactions, to minimise pipetting error. Avoid exposing fluorescently labelled probes to light.

Note: The number of reactions performed simultaneously must not exceed the maximum capacity of the thermocycler.

Table 4. Reaction mix for RAID-CRC Screen A multiplex qPCR (per reaction).

Component	Final concentration	Volume/reaction
Multiplex Master Mix 4X	1X	2.50 µL
Primer SCR_A_f1	50 nM	0.20 µL
Primer SCR_A_r1	50 nM	0.20 µL
Probe SCR_A_FAM	60 nM	0.24 µL
Primer SCR_A_f2	300 nM	0.60 µL
Primer SCR_A_r2	300 nM	0.60 µL
Probe SCR_A_HEX	100 nM	0.20 µL
Primer SCR_A_f3	300 nM	0.60 µL
Primer SCR_A_r3	300 nM	0.60 µL
Probe SCR_A_ROX	150 nM	0.30 µL
RNase-free water	-	1.96 µL
Total reaction volume		8 µL

Mix the reaction vigorously and briefly centrifuge. Dispense **8 µL** of the reaction mix into qPCR tubes compatible with the thermocycler in use.

Add **2 µL of the DNA sample** to each tube containing the reaction mix. Also, add **2 µL of the specific positive control** for Screen A to the designated control tube, and leave one tube containing only the reaction mix as the **NTC**.

Seal the tubes with optical caps, vortex for 5 seconds, and briefly spin down to ensure the reaction settles at the bottom of the tube without droplets or bubbles.

4. Multiplex qPCR Analysis – Screen B:

Add the components listed in Table 5 into a microcentrifuge tube. It is recommended to prepare a reaction mix volume of $n \times 1.1$, where n is the number of reactions, to minimise pipetting error. Avoid exposing fluorescently labelled probes to light.

Note: The number of reactions performed simultaneously must not exceed the maximum capacity of the thermocycler.

Table 5. Reaction mix for RAID-CRC Screen B multiplex qPCR (per reaction).

Component	Final concentration	Volume/reaction
Multiplex Master Mix 4X	1X	2.50 µL
Primer SCR_B_f1	300 nM	0.60 µL
Primer SCR_B_r1	300 nM	0.60 µL
Probe SCR_B_HEX	100 nM	0.20 µL
Primer SCR_B_f2	250 nM	0.50 µL
Primer SCR_B_r2	250 nM	0.50 µL
Probe SCR_B_CY5	200 nM	0.40 µL
RNase-free water	-	2.70 µL
Total reaction volume		8 µL

Mix the reaction vigorously and briefly centrifuge. Dispense **8 µL** of the reaction mix into qPCR tubes compatible with the thermocycler in use.

Add **2 µL of the DNA sample** to each tube containing the reaction mix. Also, add **2 µL of the specific positive control** for Screen B to the designated control tube, and leave one tube containing only the reaction mix as the **NTC**.

Seal the tubes with optical caps, vortex for 5 seconds, and briefly spin down to ensure the reaction settles at the bottom of the tube without droplets or bubbles.

5. Singleplex qPCR Analysis – Screen C:

Add the components listed in Table 6 into a microcentrifuge tube. It is recommended to prepare a reaction mix volume of $n \times 1.1$, where n is the number of reactions, to minimise pipetting error. Avoid exposing fluorescently labelled probes to light.

Note: The number of reactions performed simultaneously must not exceed the maximum capacity of the thermocycler.

Table 6. Reaction mix for RAID-CRC Screen C singleplex qPCR (per reaction).

Component	Final concentration	Volume/reaction
Multiplex Master Mix 4X	1X	2.50 µL
Primer SCR_C_f1	300 nM	0.60 µL
Primer SCR_C_r1	300 nM	0.60 µL
Probe SCR_C_FAM	250 nM	0.50 µL
RNase-free water	-	3.80 µL
Total reaction volume		8 µL

Mix the reaction vigorously and briefly centrifuge. Dispense **8 µL** of the reaction mix into qPCR tubes compatible with the thermocycler in use.

Add **2 µL of the DNA sample** to each tube containing the reaction mix. Also, add **2 µL of the specific positive control** for Screen C to the designated control tube, and leave one tube containing only reaction mix as the **NTC**. Seal the tubes with optical caps, vortex for 5 seconds, and briefly spin down to ensure the reaction settles at the bottom of the tube without droplets or bubbles.

6. Loading the thermocycler:

Select the appropriate channels (targets/dyes) for fluorescence acquisition during the combined hybridisation/extension step: **FAM, HEX, ROX** for multiplex **SCR_A**, **HEX, CY5** for multiplex **SCR_B** and **FAM** for singleplex **SCR_C**.

Table 7. Thermocycler protocol for RAID-CRC Screen multiplex/singleplex qPCR.

Step		Time (min:s)	Temperature (°C)
qPCR activation		01:00	95
40 cycles	Denaturation	00:15	95
	Annealing + Extension	00:30	60

7. Starting the analysis (run):

For each qPCR assay (SCR_A, SCR_B, and SCR_C), all samples and controls of the same type must be **analysed within the same run**. If sample numbers allow, **it is recommended to perform all three analyses in a single run**.

Note: If different thermocyclers are used for analyses of the same sample (e.g., SCR_A on thermocycler 1 and SCR_B on thermocycler 2), **both instruments must be the same model to ensure consistent results**.

Analysis and interpretation of the results:

1. Data processing:

Sample analysis is performed using the qPCR instrument software, following the manufacturer's instructions.

Note: Before beginning the analysis, predefined analysis parameters (baseline and threshold values) must be selected for each primer/probe system, as indicated in the Technical Specifications of the RAID-CRC Screen kit. These are provided upon kit purchase and available at the GoodGut Professional Area (<https://professionalarea.goodgut.eu/>).

2. Uploading results to the platform:

To obtain the RAID-CRC Screen diagnostic output, the results from all three qPCR analyses (multiplex and singleplex), including positive and negative controls, must be uploaded to the **GoodGut-Test™** web platform (<https://goodgut-test.eu/>) following the **User Manual**.

Data must be uploaded using Excel templates specific to each multiplex assay, which must include:

- Sample identifier
- Detection channel (dye)
- Raw Ct (Cq) value

The Excel template files can be downloaded directly from the platform.

For any technical support request or feedback, please contact support@goodgut.eu

In the event of an incident—defined as any malfunction or problem occurring with this In Vitro Diagnostic Medical Device during use or thereafter, which may have serious consequences for health—please contact the manufacturing laboratory (GoodGut S.L.U.). E-mail: vigilance@goodgut.eu and/or the competent authority of the country where the user and/or patient is established.

Symbol description:



Contains a substance hazardous to health. Refer to safety information



In vitro diagnostic medical device in accordance with Regulation (EU) 2017/746



Keep away from moisture



Temperature limits



Read the instructions for use before using the product



Number of reactions that can be performed with the kit



Date of manufacture



Reference number



Lot number



Date of expiration

RAID-CRC Screen qPCR Kit
Colorectal Cancer Screening qPCR Kit

Basic UDI-DI: 8437023437RAIDCRCKC

ANNEX 1: Compatibility of DNA extraction kits and automated equipment

The DNA extraction kits and automated extractors that may be used to obtain reliable diagnostic results with the RAID-CRC Screen assay are listed below (not included in the kit).

DNeasy PowerSoil Pro DNA Extraction Kit from QIAGEN ([manual extraction](#))

- Kit reference: 47014 (50 reactions); 47016 (250 reactions), QIAGEN.
- Before initiating DNA extraction, the sample preparation process described in the ***Pre-analytical requirements*** section of ***Requirements for the RAID-CRC Screen analysis*** must be performed. Once this preparation is complete, proceed with DNA extraction according to the manufacturer's instructions.

QIAcube from QIAGEN ([automated extraction](#))

- Use the DNeasy PowerSoil Pro DNA Extraction Kit (QIAGEN) with the QIAcube Connect automated extractor (QIAGEN).
- Before initiating DNA extraction, the sample preparation process described in the ***Pre-analytical requirements*** section of ***Requirements for the RAID-CRC Screen analysis*** must be performed. Once this preparation is complete, proceed with DNA extraction according to the manufacturer's instructions.



ANNEX 2: Compatibility of real-time PCR equipment

The RAID-CRC Screen multiplex/singleplex assays can be performed on all qPCR thermocyclers equipped with a **low-profile block**, as listed below.

AriaDx (Agilent Technologies)

- Sample analysis is performed using the software provided with the real-time PCR instrument, following the manufacturer's instructions.
- Before analysing the data, select the **predefined analysis settings** for each primer + probe set (e.g., reference settings and threshold values) according to the **Technical Specifications of the RAID-CRC Screen qPCR Kit**.
- Use tube strips and optical caps recommended by the thermocycler manufacturer.

CFX96 (BioRad)

- Sample analysis is performed using the software provided with the real-time PCR instrument, following the manufacturer's instructions.
- Specifications for analysing results with the CFX96 software:
 - Select **BR White** as the plate type.
 - Apply **fluorescence drift correction**.
- Before analysing the data, select the **predefined analysis settings** for each primer + probe system (e.g., reference configuration and threshold values) according to the *Technical Specifications of the RAID-CRC Screen qPCR Kit*.
- Use tube strips and optical caps recommended by the thermocycler manufacturer.

Note: The RAID-CRC Screen Technical Specifications for each lot are provided separately upon purchase of the kit and are also available in the Professional Area of the GoodGut website: <https://professionalarea.goodgut.eu/>.

