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## RAID-CRC Screen qPCR Kit Colorectal Cancer Screening qPCR Kit

### QUICK START PROTOCOL

The RAID-CRC Screen qPCR Kit is an *in vitro* diagnostic device (IVD) intended for professional laboratory use by qualified personnel.

#### Intended use

The RAID-CRC Screen qPCR Kit is an *in vitro* diagnostic device intended for the screening of advanced colorectal neoplasia, including precancerous lesions and colorectal cancer, in asymptomatic individuals aged 50 to 69 years (both inclusive) who have obtained a positive result in the faecal immunochemical test (FIT) at a cut-off point of 100 ng haemoglobin per mL using the collection tube from Eiken Chemical.

The RAID-CRC Screen assay detects specific bacterial markers in DNA extracted from stool samples using quantitative polymerase chain reaction (qPCR).

#### Storage and stability

The RAID-CRC Screen qPCR Kit must be transported between 2 °C and 8 °C. Upon receipt, the Master Mix 4X and positive controls must be stored between –30 °C and –15 °C in a constant-temperature freezer protected from light. The lyophilised primers and probes can be kept at room temperature until resuspension in Tris-HCl buffer (pH 8.0); once resuspended, they must be stored between –30 °C and –15 °C.

#### Safety information

- For professional use only.
- Do not use after the expiry date.
- Design a unidirectional workflow: start in the Extraction Area, then move to the Amplification and Detection Area. **Do not return samples, equipment, or reagents to previous areas.**
- Prepare reaction mixes in a DNA-free biosafety cabinet. **Avoid loading samples and positive controls in the same cabinet.**
- Follow Good Laboratory Practice: wear protective clothing, disposable powder-free gloves, safety glasses, and a face mask. Do not eat, drink, or smoke in the laboratory. Wash your hands after handling samples.
- Dispose of consumables and qPCR reagents in biohazard waste containers.
- Regularly decontaminate all working equipment, especially pipettes and work surfaces.

**Warning:** The kit contains 1,2,4-triazole. May impair fertility or cause harm to the unborn child. May cause harm to breastfed children. Seek specific instructions before use. Do not handle the kit without reading and understanding all safety information. Follow the good laboratory practices described above and ensure all personnel are informed of the associated risks. **In case of exposure, consult a physician.**

**Warning:** Do NOT add bleach or acidic solutions directly to sample waste, as toxic irritating vapors may be produced.

#### Analysis requirements

The RAID-CRC Screen qPCR Kit is optimised for stool samples obtained from asymptomatic individuals aged 50–69 years with a positive FIT result (Eiken Chemical, cut-off 100 ng Hb/mL). If a different FIT cut-off is used, please contact [support@goodgut.eu](mailto:support@goodgut.eu).

To review all pre-analytical requirements, sample handling procedures, full exclusion criteria, the pre-extraction preparation process, and the DNA extraction workflow, please refer to the Instructions for Use (IFU).

**Results are only valid when using the extraction kits and automated extractors, as well as the thermocyclers, that have been validated for this assay (see Annexes 1 and 2 of the IFU).**

#### Analysis protocol

- **qPCR procedure:** The Master Mix 4X and the pre-loaded primers/probes, as well as the parameters (temperature, number of cycles, and step times), have been optimised for optimal performance and specificity. **Note:** *Before starting, resuspend primers and probes in Tris-HCl pH 8.0, incubate for 1 h at room temperature or overnight at 4 °C, and then store between –30 °C and –15 °C.*

[EN]

**For each sample, perform two multiplex qPCRs and one singleplex qPCR: Screen A (SCR\_A), Screen B (SCR\_B), and Screen C (SCR\_C). For each qPCR analysis (SCR\_A, SCR\_B, SCR\_C), perform steps 1–5 separately.**

1. Determine the number of microtubes and caps required to perform the necessary reactions, taking into account the samples and the controls for each qPCR assay (material not included in the kit). A positive control and a no-template control (NTC) must be included in every qPCR assay. **Note:** *Each qPCR assay has its own specific positive control.*
  2. Thaw the Master Mix 4X, primers, probes, and positive controls.
  3. Add the necessary components for each qPCR analysis to a microcentrifuge tube (see the components required for each analysis in the IFUs). It is recommended to prepare a Mix volume of  $n \times 1.1$  (where  $n$  is the number of reactions) to minimise pipetting error. Minimise exposure of the fluorescently labelled probes to light. **Note:** *The number of simultaneous reactions must not exceed the capacity of the thermocycler.*
  4. Mix thoroughly, briefly centrifuge, and **dispense 8 µL into the qPCR tubes** recommended by the thermocycler manufacturer.
  5. **Add 2 µL of the DNA sample** to each reaction tube. **Add 2 µL of the specific positive control** to its designated tube and leave one tube with the reaction mix only as **NTC**.
  6. Close the tubes with the appropriate **optical caps, vortex briefly (5 seconds)**, and perform a short spin to ensure that the reaction mix settles at the bottom of the tube, avoiding the formation of droplets and/or bubbles.
  7. Insert the strips into the compatible thermocycler (see compatibility in Annex 2 of the Instructions for Use).
  8. Program the thermocycler according to the cycling conditions described in Table 7 of the IFUs. All three qPCR analyses follow the same program. **Note:** *Select fluorogenic channels during the combined annealing/extension step: FAM, HEX and ROX for Screen A; HEX and CY5 for Screen B; FAM for Screen C.*
  9. For each qPCR assay (SCR\_A, SCR\_B, or SCR\_C), all samples and controls of the same type must be analysed within the same run. In addition, it is recommended to perform the three multiplex analyses in a single run to ensure consistent experimental conditions and to facilitate interpretation of the results.
  10. Insert the qPCR tubes into the real-time thermocycler and start the run. **Note:** *If different thermocyclers are used for the analyses of the same sample (e.g., Multiplex 1 in thermocycler 1 and Multiplex 2 in thermocycler 2), both devices must be of the same model to ensure result consistency.*
- **Analysis and interpretation of results:**
    1. Sample analysis must be carried out using the software provided with the qPCR instrument, following the manufacturer's instructions for use. **Note:** *Before data analysis, select the predefined parameters (baseline and threshold) for each primer/probe set as specified in the Technical Specifications of the RAID-CRC Screen. This information is provided with the kit and available in the Professional Area of the GoodGut website (<https://professionalarea.goodgut.eu>).*
    2. To obtain the RAID-CRC Screen diagnostic result, the data generated in each multiplex qPCR assay (including positive and negative controls) must be uploaded to the **GoodGut-Test™** web platform (<https://goodgut-test.eu>) following the **User Manual**.  
The data must be uploaded using the specific Excel templates for each multiplex assay, which must include sample identifier, detection channel (dye) and raw Ct (Cq) value. The Excel templates can be downloaded directly from the platform as indicated in the User Manual.

## Further information

- **Instructions for Use (IFU) of the RAID-CRC Screen qPCR Kit:** <https://professionalarea.goodgut.eu>.  
If the website is temporarily unavailable and the IFU cannot be accessed, it can be requested by email at [support@goodgut.eu](mailto:support@goodgut.eu). The instructions will be sent within 7 calendar days from receipt of the request.  
**This quick-start protocol does not replace the Instructions for Use.** The user must consult the IFU before using the product to ensure proper handling and regulatory compliance.
- **Technical assistance:** [support@goodgut.eu](mailto:support@goodgut.eu)

Scan the QR code to access the Instructions for Use.

**GOODGUT S.L.U.**

CIF/NIF: B55206916 | GoodGut SRN: ES-MF-000000229

C/ Pic de Peguera, 11, 17003 Girona - Cataluña, España

Teléfono: +34 972 18 32 20 | E-mail: [info@goodgut.eu](mailto:info@goodgut.eu) | [www.goodgut.eu](http://www.goodgut.eu)

QG-RAID-CRC-SCR-001 v2.0

**Revision date:** Apr 2026

The information provided in this document may be subject to change due to continuous technological updates.